



# Slim Line Easy Pneumatic Crown and Bridge Remover II

**IMPROVED!**  
Smaller  
All Stainless Steel  
With Quiet Air



- Contains:
- Handpiece
  - 1 Interproximal bridge
  - 1 Single crown straight
  - 1 Single crown contra-angle
  - Open -end adjustment wrench
  - Flexi-wire set
  - Instructional/Clinical CD
  - Storage/Carrying case
  - One year limited warranty on parts and labor

- Predictable** – consistent and reliable time after time
- Safe** – Less chance of tooth fracture than any other C&B removers
- Economical** – restorations removed intact and reusable
- High Patient Acceptance** – no shocks to the patient head
- Time Saving** – significantly reduce chair time
- Easy to use** – finger tip control, no spring to reset and no drilling involved
- Maintenance Free** – no costly repairs or “down time”
- OSHA Compliance** – completely autoclavable
- No Anesthesia Required** – an atraumatic procedure
- Guarantee** – 30 Days Money Back
- Can be used without an assistant**

*“This easy remover delivers a consistent impact force. The secret of removing a crown or a bridge is the Speed; very little force is needed to break the cement seal.”*

**Dr. Carl Misch**  
Misch Implant Institute, Michigan  
Not a paid endorse

- Replacement Tips**  
Specify:
- |                           |        |
|---------------------------|--------|
| Interproximal Bridge      | Tip #1 |
| Single-Crown Straight     | Tip #2 |
| Single-Crown-Contra-Angle | Tip #3 |

According to Dr. Francis V. Panno, Professor and Division Head of the Department of Restorative & Prosthetic Sciences at New York University College of Dentistry, The Easy Pneumatic Crown and Bridge Remover will break the seal of any type of bonding cement. Elastic cements tend to retract back to the starting point when time spent recalibrating an instrument or preparing it for another “strike”. The rapid and repeated action of the Easy Pneumatic Remover prevents the cement from reseating itself.

