

# EasyScan 3-in-1 DENTAL X-RAY SCANNER

**EasyScan** is a compact electronic device that allows three different operations.

- Scan dental x-ray film for viewing on a display screen
- Scan dental x-ray film directly to a computer to archive, email, etc.
- Connect intra-oral camera or other analog video devices directly to a computer

**EasyScan** provides on-touch scanning of standard 31x41mm (1.2"x1.6") dental x-ray film enabling it to be magnified and viewed on a display screen. At the same time, you can choose to save the image to a computer for archival, email or other purposes.

**EasyScan** can also be used to connect an analog intra-oral camera to a computer for direct snapshot storage.

**EasyScan** is an indispensable device for every dental office and clinic.