

POWER Light

Best choice for the best vision



POWER Light is a new dental light of Admetec Solutions Ltd. LED light developed based on a newest technologies and the last world achievements.

The Results:

- Highly powerful light beam (> 70,000lx)
- Very low weight (12g)
- High light concentration

POWER Light creates homogeneous square light beam. Like the previous models of Admetec lights it's powered from PowerPack - intelligent rechargeable battery pack.



Regular light



Power Light

POWER Light produces in two models: "headband" and "loupes". It has attachments to most of known loupes brands: Orascoptic, Heine, SurgiTel, Zeiss and other.



Technical specifications

Lightness	- >70,000lx
Light source weight	- 12G
Light color temp	- 6000K
Battery working time	- 5hours
Recharge time	- 2hours
LED life time	- 50,000h