

OrBo

Extreme Oral Cavity Illumination

The OrBo system is tiny, yet powerful light source for the illumination of an oral cavity while treating you patient.

The OrBo is a mini U-shaped clip equipped with 2 high brightness LEDs. The clip which comes in two sizes one for suction tub and one for the saliva ejector are autoclavable.

Once the system is on the OrBo will get you double the amount of light you usually get in the oral cavity.

The long tiny power cord is the “Brain” behind the light; the cord consists of a unique electronic design, stabilizing the light out put. The cord easily attached to the ejector hosing with the provided Velcro stripes thus not interfering with work.

The OrBo is a hands free powerful light source, easy to operate that will ease eye stress.